

General Description	This centering device should allows to quickly carry out hardness test of O rings with cross section diameter between 1.5mm to 11mm.
Oring dimension to be tested	Device should be capable of holding Oring with outer diameter of 5mm to 200 mm and Oring cut pieces of minimum 10mm length.
Fixing mechanism	The device has to be applied to the standard plate of the instrument through magnetic fixing which allows rapid installation and removal The adjustment wheel allows you to set the distance between the cylinders according to the cross section diameter of the Oring. The fixing clip for the correct positioning of curved pieces. The integrated extension plate should allow holding large pieces of O ring up to 200mm ID.